


EYEZEN[®] LENSES CUTOUT CHARTS

EYEZEN® START™ & EYEZEN®+ 1, 2, 3 AND 4 CUTOUT CHART

Eyezen+: Fitting Height required with minimum of 15mm
Eyezen Start: Fitting Height required with no minimum value



Please verify correct size with a ruler

¹Symbol location for Lens Design - ◇

Symbol location for Lens Material - M

| LENS DESIGN LENS ENGRAVING SYMBOLS | |
|---------------------------------------|--------|
| PRODUCT | SYMBOL |
| Eyezen+ 1 (0.40 D) | ⊕ |
| Eyezen+ 2 (0.60 D) | ⊕ |
| Eyezen+ 3 (0.85 D) | ⊕ |
| Eyezen+ 4 (1.10 D) | ⊕ |
| Eyezen Start | ⊕ |

| LENS MATERIAL LENS ENGRAVING SYMBOLS | |
|---|--------|
| MATERIAL | SYMBOL |
| Plastic 1.50 | |
| TREXA® (Trivex®) 1.53 | ✳ |
| Airwear® Polycarbonate 1.59 | P |
| Thin & Lite® 1.60 | 6 |
| Thin & Lite® 1.67 | Π |
| Thin & Lite® 1.74 | 74 |


PRINTING INSTRUCTIONS

To print this measuring chart at the correct size from Acrobat Reader:

1. Go to **File**, then select **Print**
2. Under **Page Sizing & Handling** options, select **Actual Size**
3. Click on **Print** button to print the chart.

Note: Chart measurements will be incorrect if "Fit" or other page-scaling options are selected. It is recommended to verify that the printed chart is the correct size with a ruler.

EYEZEN® + 0 CUTOUT CHART



Please verify correct size with a ruler

¹Symbol location for Lens Design -

Symbol location for Lens Material - **M**

| LENS DESIGN LENS ENGRAVING SYMBOLS | |
|---------------------------------------|--------|
| PRODUCT | SYMBOL |
| Eyezen+ 0 (0.00 D) | |

| LENS MATERIAL LENS ENGRAVING SYMBOLS | |
|---|--------|
| MATERIAL | SYMBOL |
| Plastic 1.50 | |

TREXA®
(Trivex®) 1.53

Airwear®
Polycarbonate 1.59

Thin & Lite® 1.60

Thin & Lite® 1.67

Thin & Lite® 1.74

PRINTING INSTRUCTIONS

To print this measuring chart at the correct size from Acrobat Reader:

1. Go to **File**, then select **Print**
2. Under **Page Sizing & Handling** options, select **Actual Size**
3. Click on **Print** button to print the chart.

Note: Chart measurements will be incorrect if "Fit" or other page-scaling options are selected. It is recommended to verify that the printed chart is the correct size with a ruler.

BLOCKING AND FINISHING

Eyezen® enhanced single vision lenses are an everyday single vision lens; however, Eyezen® Start™, Eyezen+ 1, 2, 3 and 4 designs should be blocked and finished similarly to a progressive lens:

- FRP (Fitting Cross 4mm above the plane of the micro circles) should be used for **blocking for finishing**.

For Eyezen+ 0 design, fit as a standard single vision lens.

Eyezen static cling verification decals (12593_PRO_ZEN) are available from Essilor.

VERIFICATION

Instructions to verify power in
Eyezen Start and *Eyezen+ 1, 2, 3* and *4* lenses

- 1 Verify the prescribed power at the distance verification circle
- 2 If prism was prescribed, then verify prism power at PRP
- 3 Confirm engraving to ensure the correct *Eyezen* lens was provided

Note: No other checks are necessary to confirm accommodative relief on *Eyezen+ 1, 2, 3* and *4* lenses

Instructions to verify power in
Eyezen+ 0 lenses

- 1 Verify the prescribed power at the fitting cross
- 2 If prism was prescribed, then verify prism power at the fitting cross
- 3 Confirm engraving to make sure an *Eyezen+ 0* lens was provided

Patients will not receive the full benefit of the accommodative relief in *Eyezen+ 1, 2, 3* or *4* if the fitting height is less than 15mm. Fitting height is also required on *Eyezen Start* with no minimum value. For *Eyezen+ 0* design, including a fitting height is recommended as *Eyezen+ 0* design is a digital aspheric single vision lens.