

## Troubleshooting Fitting Issues


Distance vision is

slightly blurry.

**SYMPTOM** SOLUTION Verify fitting height and PD Patient has narrow measurements reading area. Verify ADD power Add pantoscopic tilt and decrease vertex distance Adjust frames to decrease Peripheral vision blurs vertex distance and to and moves. increase facial wrap Lenses are fit too low: Patient lifts head or glasses to read. Adjust frames to sit higher on patient's face Adjust nose pads closer together Decrease pantoscopic tilt If necessary, refit lenses Lenses are fit too high: Patient lowers head Adjust frames to sit lower or glasses to read on patient's face at a distance. Lower frames by widening nose pads Increase pantoscopic tilt If necessary, refit lenses PD is off or lenses are mounted Patient moves reading incorrectly: material off to side for Use lens cling to verify better focus. lenses are on iris Verify monocular PD measurements Have lenses remade with correct PD measurements

Increase pantoscopic tilt

Decrease fitting height



FITTING & **DISPENSING GUIDE** VARILUX® Progressive Lenses **essilor VariluX**<sub>®</sub>

## **Fitting**

## Dispensing

## **Understanding**

### **ADJUST THE FRAME**

- Adjust the frame on the patient for maximum comfort and accuracy before taking any measurements
- Set the vertex distance between 12 and 14mm
- Set the pantoscopic tilt angle between 6° and 8°
- Frame should have positive facial wrap



### **MEASURE PUPILLARY DISTANCE (PD)**

- Always take monocular PD to ensure exact centering of the eye behind the lens
- Varilux<sup>®</sup> lenses should be fitted using distance monocular PD




### **MEASURE FITTING HEIGHT**

- Situate yourself directly in front of patient's eyes
- Mark each lens at center pupil using a felt-tip pen
- Draw a horizontal line on each lens and double-check to make sure that the lines are crossing the center of each pupil
- Measure the fitting height from the deepest point of the lens to the horizontal line
- · For zyl frames, measure to the lowest point of the bevel

# 18mm


### **CHECK LENS CUT-OUT**

- Mark the patient's fitting height and distance PD on the sample lens, creating a cross
- Place the lens cross over the layout chart cross to verify that the lens will fit within the lens diameter circle
- If the lens does not fit, choose another, more suitable frame that will accommodate the lens



### **CONFIRM MEASUREMENTS AND PRESCRIPTION**


- Using the centering chart, center the frame over the inverted "V"
- · Confirm the monocular PD and fitting height
- · Confirm distance Rx with lensometer\*
- Confirm prism at prism reference point (PRP)



\*The abbreviated ADD power is engraved on the temporal side of all Varilux lenses and should be used to verify ADD power. Abbreviated ADD power is the first two digits of the actual ADD power to the nearest 0.25 diopters. For example, 22 is equal to 2.25.

### **RECREATE FITTING CROSS (IF NEEDED)**

- Mark the lens micro-circles with a felt tip pen
- Lay the glasses over the cut-out chart and align the micro-circles to the micro-circles on the cut-out chart
- Mark the fitting cross with a felt tip pen




### **VERIFY COMPENSATED Rx WITH A DUAL PACK SLIP**


- Verify the original prescription matches what was ordered
- Apply ANSI tolerances to compensated values

	Rx Ordered	Sphere	Cylinder	Axis	ADD
R	Vx Phys W3+	+0.00	-0.50	180	3.00
L	Vx Phys W3+	+0.00	-0.50	180	2.25
	Rx in Lensometer	Sphere	Cylinder	Axis	ADD
R	Rx in Lensometer Vx Phys W3+	Sphere 0.04	Cylinder -0.46	Axis 180	ADD 2.84

### **CONFIRM FIT ON PATIENT**

- With lenses marked or using decals, verify that the fitting of cross is at center pupil
- Adjust the frame to raise or lower the fit, if necessary





The following markings are used to verify the properties of a progressive addition lens (PAL):

- 1 Fitting Reference Point (FRP) Also called the fitting cross. Designates the point on the lens that should be placed along the optical axis of the patient (center of pupil). Once removed, this marking may be recreated using a layout chart.
- **2 Prism Reference Point (PRP)** Located 4mm below the FRP, this is the optical center of a PAL, and it's used to check prismatic properties. Prescribed or thinning prism should be verified at this point.
- **3 Two Engraved Micro-Circles** Situated 17mm to each side of the PRP, these circles are used to verify the axis alignment and can be used with a layout chart to recreate the FRP marking.
- **4 Design Identifier** An engraving located under the nasal micro–circle that's unique to each PAL design.
- **5** Material Identifier Appears to the right of the design identification and denotes the material.
- **6 ADD Engraving** Appears under the temporal micro–circle and denotes the ADD power. Verification of ADD power should be made with the engraving, not with the near verification circle.
- **7 Distance Verification Circle** Located above the FRP, this area of the lens may be used to verify prescribed distance power but cannot be used to verify prism.
- **8** Near Verification Circle Located 13–18mm below the FRP, the near verification circle is no longer used to verify total near power. Modern PALs have varied insets and progression lengths that may place the near area outside of the circle. True ADD power cannot be read by a lensometer, since the worn position of the lens will vary from the lensometer mounting.