

A man with short brown hair and glasses, wearing a dark blue suit jacket over a light blue collared shirt, is smiling and looking slightly to the right. The background is a blurred cityscape with a prominent white wind turbine tower with yellow bands. The overall tone is professional and positive.

**VARILUX®**

**LENSES AVAILABILITY**  
for Eyecare Professionals

Effective **AUGUST 11, 2020**



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX® X DESIGN		LENS ID	VARILUX® X FIT											LENS ID	VARILUX® X 4D			LENS ID
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER			
PLASTIC	• Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.53	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00			
THIN & LITE®	• Clear	1.67											-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.67	✓SM	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00			
	• Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™***	1.74	✓	✓									-14.00 to +13.00	0.75 to 4.00D	-6.00			

MFH=14mm

VARILUX® S DESIGN		LENS ID	VARILUX® S FIT											LENS ID	VARILUX® S 4D			LENS ID
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER			
PLASTIC	• Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Index Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.53	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.59	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00			
THIN & LITE®	Clear	1.67											-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive®	1.67	✓	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00			
	Transitions® Signature® GEN 8™***	1.74	✓	✓									-14.00 to +13.00	0.75 to 4.00D	-6.00			

MFH=14mm

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

• **G:** Gray, **B:** Brown, **GG:** Graphite Green, **Colors:** Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, **G1:** Gray 1, **B1:** Brown 1, **G-G:** Gray-Green, **C:** Copper, **R:** Ruby, **ESC:** Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX® PHYSIO® W3+		TEMPORAL	NASAL	VARILUX® PHYSIO® W3+ FIT										TEMPORAL	NASAL	VARILUX® PHYSIO® W3+ EYECODE™				TEMPORAL	NASAL
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER					
PLASTIC	• Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Essilor Polarized	1.50	✓	✓			✓					☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Xperio UV™ Polarized	1.50	✓	✓			✓					☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® XTRActive®	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00					
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Essilor Polarized	1.59	✓	✓			✓					☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Xperio UV™ Polarized	1.59	✓	✓			✓					☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00					
	Transitions® XTRActive®	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00					
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Transitions® Signature® GEN 8™	1.60	✓	✓	✓								-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00					
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00					
Transitions® Signature® GEN 8™	1.74	✓										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00						

MFH=14mm

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, COLORS: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

▼ Availability date to be announced. Current availability Transitions® Signature® VII Only.

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX® COMFORT MAX		LENS ID										VARILUX® COMFORT MAX FIT		LENS ID									
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER							
PLASTIC	• Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Transitions® XTRActive®	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
	Transitions® XTRActive®	1.53	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00							
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00							
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00							
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00							
	Transitions® XTRActive®	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00							
THIN & LITE®	• Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00							
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00							
	• Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00							
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00							
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00							
	Transitions® XTRActive®	1.67	✓SM										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00							

**MFH=14mm**

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX COMFORT® W2+		TEMPORAL NASAL												VARILUX COMFORT® W2+ FIT		TEMPORAL NASAL												VARILUX COMFORT® W2+ EYECODE™		TEMPORAL NASAL	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER															
PLASTIC	• Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Transitions® XTRActive®	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00															
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00															
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00															
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00															
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00															
	Transitions® XTRActive®	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00															
THIN & LITE®	• Clear	1.60				✓							-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00															
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00															
	• Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00															
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00															
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00															
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00															
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00															

**MFH=14mm**

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

\*\**THIN & LITE* 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

## VARILUX PHYSIO DRx™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	• Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓	✓							-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.50	✓										-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

▼ Availability date to be announced. Current availability Transitions® Signature® VII Only.

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

## VARILUX PHYSIO SHORT DRx™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® Vantage®	1.50	✓										-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.00D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00
Transitions® Signature® GEN 8™	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00	

MFH=14mm

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*\**THIN & LITE* 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

▼ Availability date to be announced. Current availability Transitions® Signature® VII Only.

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX COMFORT DRx™		LENS ID										POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER		
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R					ESC	
PLASTIC	Clear	1.50												-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00
	Essilor Polarized	1.50	✓	✓				✓	✓	✓		☼		-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓				✓	✓	✓		☼		-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓									-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
	Transitions® Vantage®	1.50	✓											-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX)	• Clear	1.53												-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59												-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓				✓	✓	✓	✓	☼		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓				✓	✓	✓	✓	☼		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓									-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00
	Transitions® Vantage®	1.59	✓											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00
THIN & LITE®	Clear	1.60												-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓										-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	Clear	1.67												-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.67	✓	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	Clear**	1.74												-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00
Transitions® Signature® GEN 8™	1.74	✓											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	

MFH=17mm

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

▼ Availability date to be announced. Current availability Transitions® Signature® VII Only.

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)





# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

## VARILUX COMFORT SHORT DRX™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® Vantage®	1.50	✓										-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00
TREXA® (TRIVEX)	• Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00

**MFH=14mm**

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ **Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)**



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX PHYSIO ENHANCED™		VARILUX PHYSIO ENHANCED FIT™										VARILUX PHYSIO ENHANCED EYECODE™				
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX)	Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	✓										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear	1.67											-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear*	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

**NOTES:** Varilux Physio Enhanced Eyecode lenses are exclusive to the Visiooffice® System.

\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia™, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

At fitting height less than 17mm, ADD power is up to +3.00.

G: Gray, B: Brown, GG: Graphite Green, COLORS: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

VARILUX® PHYSIO®		LENS ID H → e											POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC						
PLASTIC	• Clear	1.50												-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓										-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓										-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	Clear	1.53												-10.00 to +6.00	78/83	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓									-10.00 to +6.00	76/81	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59												-10.00 to +6.00	80/85	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓										-10.00 to +6.00	75/80	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓										-10.00 to +6.00	75/80	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	✓											-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	Clear	1.60												-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	• Clear	1.67												-12.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear*	1.74												-18.00 to +9.00	68 78/84 73/79	1.75, 3.25, 5.25, 6.50 8.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

NOTE: \*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia™, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors  
 ✓ **Currently Available**

VARILUX PHYSIO SHORT™		LENS ID → e S											POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC						
PLASTIC	• Clear	1.50												-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
TREXA® (TRIVEX®)	Clear	1.53												-10.00 to +6.00	81/86	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59												-10.00 to +6.00	80/85	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.00D	-4.00
THIN & LITE®	Clear	1.60												-12.00 to +8.00		1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.00D	-4.00
	• Clear	1.67												-12.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.00D	-4.00

MFH=14mm

• Available in Essential Blue Series™ lenses.

✓ **Currently Available**



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

## VARILUX COMFORT ENHANCED™



## VARILUX COMFORT ENHANCED FIT™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.50	✓										-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	✓										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00

MFH=14mm

NOTE: G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors  
 ✓ Currently Available



# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **AUGUST 11, 2020**

## VARILUX COMFORT®



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	• Clear	1.50											-10.00 to +7.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
	Transitions® XTRActive®	1.50	✓										-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Transitions® XTRActive®	1.59	✓										-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	73/78	1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.50D	-6.00
	• Clear	1.67											-14.00 to +9.00	78/83	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.50D	-6.00

**MFH=17mm** • Available in Essential Blue Series™ lenses.

## VARILUX COMFORT SHORT™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	• Clear	1.50											-10.00 to +7.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
	Transitions® XTRActive®	1.50	✓										-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Transitions® XTRActive®	1.59	✓										-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.00D	-4.00
	• Clear	1.67											-14.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.00D	-4.00

**MFH=14mm**


**NOTE:** At fitting height less than 17mm, ADD power is up to +3.00.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

■ Available in Transitions Signature GEN 8 Q1 2020

✓ Currently Available

# Varilux® Digital Progressive Lenses



ESSILORUSA.COM