

CRIZAL® RX AND STOCK LENS

AVAILABILITY CHART

CRIZAL®
#1

PERFORMING NO-GLARE BRAND
ACCORDING TO INDEPENDENT
EYECARE PROFESSIONALS¹

NEW

Crizal® SapphireTM HR

NO-GLARE LENSES

THE BEST
ANTI-REFLECTIVE
COATING OVERALL²

¹ Survey conducted in 2020 by an independent third party, sponsored by Essilor.
383 Independent Eyecare Professionals were surveyed.

² Within the US market, based on the results in several studies in areas of anti-reflective efficiency, ease of cleaning, smudge and scratch resistance, and UV protection
(Internal and External testing funded by Essilor, 2020).

Frames: Ray-Ban 3447-N ROUND METAL



OCTOBER 2022

TABLE OF CONTENTS

Crizal® Rx Lens Availability

Varilux® progressive lenses	3
DEFINITY® and Essilor PALs progressive lenses	4
Eyezen® Single Vision, Essilor Single Vision, Bifocals, Trifocals, Mid Distance/Visual Fatigue and Crizal® for Kids lenses	5
Shamir Insight Inc. lenses	6
Shore lenses	7
Signet Armorlite lenses	8
KBco lenses	9
Younger Optics, X-Cel lenses and Specialty Items	10
<i>Crizal</i> Finished Single Vision stock lenses	11, 12, 13
<i>Crizal</i> lens benefits	14
<i>Crizal</i> product guide	15

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

PLASTIC 1.50		TREXA®				AIRWEAR®				THIN & LITE®																		
		Transitions® ⁴		Transitions® ⁴		Transitions® ⁴		Transitions® ⁴		1.60		1.67 ²		1.74 ³														
Clear ¹	Essilor Polarized ^{*+▲}	Signature GEN 8™	XTRACTIVE® New Generation	XTRACTIVE® Polarized™	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRACTIVE New Generation	XTRACTIVE Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRACTIVE New Generation	XTRACTIVE Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRACTIVE New Generation	XTRACTIVE Polarized	Clear ³	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRACTIVE New Generation	XTRACTIVE Polarized	Clear ³	Signature GEN 8	XTRACTIVE New Generation	XTRACTIVE Polarized
PROGRESSIVE LENSES																												
Varilux® X 4D	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® X Fit	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® X Design	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio® W3+ eyecode™	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio® W3+	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio® W3+ FIT	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort Max	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort Max FIT	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort® W2+ eyecode™	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort® W2+	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort® W2+ FIT	■	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio Enhanced eyecode™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio Enhanced™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio Enhanced FIT™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® DRx™	■	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Short DRx™	■	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort Enhanced™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort Enhanced FIT™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort DRx™	■	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Liberty™ 3.0	■	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio®	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Physio Short™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort®	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
Varilux® Comfort Short™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

■ Available in Essential Blue Series® lenses.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

		PLASTIC 1.50				TREXA®				AIRWEAR®				THIN & LITE®																												
DEFINITY	ESSILOR PALS	Clear ¹	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Transitions ^{®4}		Clear	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Clear	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Clear	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Clear	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Clear	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™			
		DEFINITY®3 PLUS	DEFINITY®3	Signature® GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive® New Generation	XTRActive® Polarized™									
DEFINITY®3 PLUS		●		●				●		●				●		●		●				●			●		●		●		●		●		●							
DEFINITY®3		●			●			●			●				●		●		●			●			●		●		●		●		●		●							
DEFINITY®		●	●	●	●			●		●				●		●		●			●			●		●		●		●		●		●		●						
DEFINITY SHORT™		●	●	●	●			●		●				●		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal® Fit™		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal® Advanced Fit™		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal®		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal Short™		●	●	●	●			■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal Advanced™		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal Advanced Wrap™		●	●	●	●	●	●			●						●		●			●			●		●		●		●		●		●		●						
Essilor Ideal® Sport		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal® Driving		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Essilor Ideal® Computer™		●	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●						
Accolade®		●		●	●	●	●		●					●		●		●			●			●		●		●		●		●		●		●						
Accolade Freedom™		●		●	●	●		●		●				●		●		●			●			●		●		●		●		●		●		●						
SmallFit®		●		●						●				●		●		●			●			●		●		●		●		●		●		●		●				
SmallFit® Digital		●	●	●	●	●		■		●				■		●		●			●			●		●		●		●		●		●		●		●				
Ovation®		●		●						●				●		●		●			●			●		●		●		●		●		●		●		●				
Ovation® Digital		●	●	●	●	●		■		●				■		●		●			●			●		●		●		●		●		●		●		●				
Essilor Natural®		●		●						●				●		●		●			●			●		●		●		●		●		●		●		●				
Essilor Natural® Digital		●	●	●	●	●		■		●				■		●		●			●			●		●		●		●		●		●		●		●				
Adaptar®		●								●				●				●																								
Adaptar® Digital		●		●	●	●		■		●				■		●		●			●			●		●		●		●		●		●		●		●				
Adaptar® Short Digital		●		●	●	●		■		●				■		●		●			●			●		●		●		●		●		●		●		●				
Brillenz™		■	●	●	●	●	●	■		●				■		●		●			●			●		●		●		●		●		●		●		●				

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

■ Available in Essential Blue Series™ lenses.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

		PLASTIC 1.50		TREXA®			AIRWEAR®			THIN & LITE®								
SINGLE VISION	BIFOCALS	Clear ¹	Essilor Polarized ^{*+▲}	Transitions ^{®4}		Transitions ^{®4}		Clear	Essilor Polarized ^{*+▲}	Transitions ^{®4}		1.60		1.67 ²		1.74 ³		
		Signature GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation
OTHER LENSES																		
Eyezen® Start Enhanced Single Vision ⁶		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Eyezen®+ Enhanced Single Vision ⁶		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Eyezen® Kids Enhanced Single Vision																		
Crizal Kids UV™ †																		
Crizal® Prevencia® Kids †																		
Single Vision		●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●
Single Vision Aspheric		●			■		●		■	●	●	●	●	●	●	●	●	●
Essilor 360°™ Single Vision					●		●	●	●	●	●	●	●	●	●	●	●	●
Essilor eyecode™ Single Vision					●		●	●	●	●	●	●	●	●	●	●	●	●
Essilor FIT™ Single Vision					●		●	●	●	●	●	●	●	●	●	●	●	●
Essilor Azio™ Single Vision					●		●	●	●	●	●	●	●	●	●	●	●	●
Essilor Stylistic® Single Vision					●		●	●	●	●	●	●	●	●	●	●	●	●
Essilor Sport Solutions™ ⁵										●								
FT28		●	●	●	●	◆		●		●	●	●	●	●	●	●	●	
FT28 Aspheric										●								
FT28 Essilor Bifocal AB™		●		●						●	●	●			●			
FT35		●		●			●			●	●	●			●		●	
FT45		●																
ROUND 22		●																
7 x 28		●	●	●						●		●			●	●		
8 x 35		●	●	●			●								●	●		
Essilor Anti-Fatigue™		●								●					●			
Essilor Computer™										●								
Interview®		●																

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

5 Essilor Sport Solutions has 6 different lens colors available in tint and polarized.

6 Eyezen, Eyezen Start, and Eyezen Kids lenses come systematic with Essential Blue Series included.

■ Available in Essential Blue Series® lenses.

* Crizal Sapphire HR available on SFSV Airwear® Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

† Crizal Kids UV and Crizal Prevencia Kids lenses are bundled products.

They come as a Single Vision polycarbonate lens, coated with the respective Crizal No-Glare lens treatment.

◆ Available in Transitions XTRActive Generation 1

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

PLASTIC 1.50				TRIVEX®				AIRWEAR®				THIN & LITE®						1.74 ³											
		Transitions ⁴				Transitions ⁴				Transitions ⁴			1.60		Transitions ⁴		1.67 ²		Transitions ⁴		1.74 ³								
Clear ¹	Essilor Polarized ^{★+▲}	Signature GEN 8™	XTRAactive®	New Generation	XTRAactive® Polarized™	Clear	Essilor Polarized ^{★+▲}	Signature GEN 8	XTRAactive®	New Generation	XTRAactive® Polarized	Clear	Essilor Polarized ^{★+▲}	Signature GEN 8	XTRAactive®	New Generation	XTRAactive® Polarized	Clear	Essilor Polarized ^{★+▲}	Signature GEN 8	XTRAactive®	New Generation	XTRAactive® Polarized	Clear ³	Signature GEN 8	XTRAactive®	New Generation	XTRAactive® Polarized	Transitions ⁴
SHAMIR	Shamir Single Vision	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph III® Fixed 11, 13, 15, 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph III Variable	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Attitude III Fashion 15 mm, 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Attitude III® Sport 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph III Single Vision	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Attitude III Single Vision	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II® Fixed 11, 13, 15, 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II+ Fixed 11, 13, 15, 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II Variable, II+ Variable	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph Intelligence™ 11,13,15,18mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph Intelligence™ Variable	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II Attitude 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II Attitude Short 15 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II Office	●																											
	Autograph II Single Vision	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Autograph II Single Vision Attitude™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Element™ 19 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Element Short 16 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Spectrum™ Single Vision	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Spectrum, Spectrum+ 14, 16, 18 mm	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Relax™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	FirstPAL™	●																											
	Golf™							●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Golf - Single Vision™							●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	InTouch™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	WorkSpace™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Computer™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Shamir Duo™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Shamir Genesis HD™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	

1 E-SPF® 35 with Crizal Sapphire® 360° UV, Crizal® Prevencia®, and Crizal Avancé UV™ / E-SPF® 25 with all other Crizal® products made with any lens material other than clear 1.50 plastic.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite® lenses come systematic with your choice of Crizal® (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

■ Available in Essential Blue Series® lenses.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia® is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

		PLASTIC 1.50			TRIVEX®			AIRWEAR®			THIN & LITE®																					
		Clear ¹	Essilor Polarized ^{*+▲}	Activations™	Transitions ³			Signature GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	1.74
		Clear ¹	Essilor Polarized ^{*+▲}	Activations™	Transitions ³			Signature GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	Essilor Polarized ^{*+▲}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ²	1.74
OTHER LENSES																																
SHORE	Shore Single Vision Spherical	●	●	●							●		●	●	●		●				●	●	●									
	Shore Single Vision Aspheric										●			●							●											
	Shoreview® ES	●		●	●						●		●	●	●		●				●											
	Shoreview Mini ES	●		●	●						●		●	●	●		●				●											
	Shoreview Digital	●	●		●						●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	Shoreview Digital Mini	●	●		●						●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	Shoreview Digital Advanced	●	●		●	●	●				●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	Shoreview Digital Advanced Mini	●	●		●	●	●	●			●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	Shoreview Digital Advanced Fit	●	●		●	●	●	●			●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	Shoreview Digital Advanced Mini Fit	●	●		●	●	●	●			●		●	●	●		●				●	●	●	●	●	●	●	●	●			
	FT28			●																												

1 E-SPF® 35 with Crizal Sapphire® 360° UV, Crizal® Prevencia®, and Crizal Avancé UV™ / E-SPF® 25 with all other Crizal® products made with any lens material other than clear 1.50 plastic.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

PLASTIC 1.50			TRIVEX®			AIRWEAR®			THIN & LITE®			1.74 ³ Younger Drivewear™ 1.50 / Polycarbonate															
Clear ¹	Essilor Polarized*▲	PhotoView™	Transitions ⁴	Signature® GEN 8™	XTRAactive® New Generation	Clear	Essilor Polarized*▲	Transitions ⁴	Signature GEN 8	XTRAactive New Generation	Clear	Essilor Polarized*▲	Transitions ⁴	1.60	Transitions ⁴	Signature GEN 8	XTRAactive New Generation	Clear	Essilor Polarized*▲	Transitions ⁴	Signature GEN 8	XTRAactive New Generation	Clear	Essilor Polarized*▲	Transitions ⁴		
KODAK Concise®	•																										
KODAK Easy®	•	•																									
KODAK Precise®	•		•	•																							
KODAK Precise® Short	•		•	•																							
KODAK Precise® PB	•	•	•	•	•	•	•																				
KODAK Precise® Short PB	•	•	•	•	•	•	•																				
KODAK Precise® Plus	•	•	•	•			•																				
KODAK Precise® Short Plus	•	•	•	•			•																				
KODAK Unique™	•	•	•	•	•	•	•																				
KODAK Unique HD™	•	•	•	•	•	•	•																				
KODAK Unique DRO®	•	•	•	•	•	•	•																				
KODAK Unique DRO HD™	•	•	•	•	•	•	•																				
KODAK SoftWear®	•																										
KODAK PowerUp®	•		•																								
KODAK Digital Single Vision	•	•	•	•	•	•	•																				
DirecTek™	•	•	•	•	•	•	•																				
DirecTek™ Short	•	•	•	•	•	•	•																				
Navigator®	•		•															•									
Navigator Short®	•																										
Navigator® FBS	•	•	•	•	•	•	•																				
Navigator Short® FBS	•	•	•	•	•	•	•																				

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

* Crizal Sapphire HR available on SFSV Airwear Colors.

▲

Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

		PLASTIC 1.50			TRIVEX®			AIRWEAR®			THIN & LITE®					
KBco	Other Lenses	Clear ¹	Essilor Polarized ^{*+▲}	Transitions ^{®3}	Clear	Essilor Polarized ^{*+▲}	Transitions ³	Clear	Essilor Polarized ^{*+▲}	Transitions ³	Clear	Essilor Polarized ^{*+▲}	1.60	1.67 ²		
		Signature GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Signature GEN 8	Transitions ³	Transitions ³
OTHER LENSES																
Single Vision		●			●			●			●			●		
Single Vision Aspheric														●		
EOS™ PAL		●														
iRx® PRO		●												●		
iRx® SHORT																
Xplore™		●						●								
FT28		●						●					●	●		
FT35		●						●					●	●		
7 x 28		●											●	●		
8 x 35		●											●	●		

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV™ available on polarized and tinted lenses.

RX LENS AVAILABILITY

for Eyecare Professionals

Effective OCTOBER 2022



Crizal®

NO-GLARE LENSES

		PLASTIC 1.50			TRIVEX®			AIRWEAR®			THIN & LITE®						THIN & LITE®						THIN & LITE®						Younger NuPolar® ⁵		
		Clear ¹	Essilor Polarized ^{*▲}	Transitions ^{®4}	Signature GEN 8™	XTRActive® New Generation	XTRActive® Polarized™	Younger Drivewear™	Clear ¹	Essilor Polarized ^{*▲}	Transitions ^{®4}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Younger Drivewear	Clear ¹	Essilor Polarized ^{*▲}	Transitions ^{®4}	Signature GEN 8	XTRActive New Generation	XTRActive Polarized	Clear ¹	Essilor Polarized ^{*▲}	Transitions ^{®4}	Signature GEN 8	XTRActive New Generation	Younger NuPolar ^{®5}				
YOUNGER OPTICS	OTHER LENSES	Younger Single Vision Spherical	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Younger FT28	●		●		●		●	●	●	●					●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	
	Younger FT35	●		●												●													●		
	Younger FT45	●																													
	Younger 7 x 28	●		●																										●	
	Younger 8 x 35	●		●																											
	Younger Image®	●		●		●		●	●	●	●					●				●	●	●	●	●	●	●	●	●	●		
	Younger Camber Lens Blank	●		●					●	●						●			●	●	●	●	●	●	●	●	●	●	●	●	
	X-CEL SFSV	●							●																						
X-CEL	X-CEL SFSV Aspheric	●							●																						
	X-CEL D-25	●																													
	X-CEL D-28	●		●	●				●	●	●																				
	X-CEL D-35	●		●					●	●	●																				
	X-CEL 7X25	●																													
	X-CEL 7X28	●		●	●				●	●	●																				
	X-CEL 7X35	●																													
	X-CEL 8X35	●		●					●	●	●																				
	X-CEL RD-22	●																													
	X-CEL RD-24	●		●					●	●	●																				
	X-CEL Blended Rd	●																													
	X-CEL Freedom ID Pal								●		●	●																			
SPECIALTY ITEMS (on UNCOATED LENSES ONLY)																															
Slaboff, High Adds, Bi-Concave, etc...																															

1 E-SPF® 35 with Crizal® Sapphire™ HR, Crizal® Prevencia®, Crizal® Rock™, Crizal® Easy Pro, and Crizal Easy®.

2 Essilor Sun Portfolio products are available in various colors, depending on material.

3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal (any Crizal coating starting with Crizal Easy).

4 Transitions products are available in Gray, Brown, Graphite Green, and more. For the complete list of Transitions color availability please refer to the Transitions availability chart.

5 Crizal SunShield™ is available with Younger NuPolar lenses in the following materials: 1.50, Trivex, 1.60, and 1.67.

* Crizal Sapphire HR available on SFSV Airwear Colors.

+ Crizal Prevencia is not available on polarized lenses.

▲ ONLY Crizal SunShield UV available on polarized and tinted lenses.



RX LENS AVAILABILITY
for Eyecare Professionals

Effective OCTOBER 2022



FINISHED Single Vision

CRIZAL® SAPPHIRE™ HR			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite® 1.74 (Aspheric)	1.74	-3.00 to -12.00 up to -2.00 cyl	70/75 mm
<i>Thin & Lite</i> 1.67 (Aspheric)	1.67	+6.00 to -12.00 up to -2.00 cyl	65/70/75/80 mm
Airwear® (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
<i>Airwear Essential Blue Series®</i> (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
<i>Airwear Transitions Signature GEN 8</i> Gray (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
Plastic (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
<i>Plastic Transitions Signature GEN 8</i> Gray (Spherical)	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70 mm

CRIZAL® PREVENCIA®			
MATERIAL	INDEX	POWER RANGE	SIZE
<i>Airwear</i> (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
Plastic (Spherical)	1.50	+6.00 to -10.00 up to -2.00 cyl	70 mm

CRIZAL® ROCK™			
MATERIAL	INDEX	POWER RANGE	SIZE
<i>Thin & Lite</i> 1.67 (Aspheric)	1.67	+6.00 to -12.00 up to -2.00 cyl	65/70/75/80 mm
<i>Airwear</i> (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
<i>Airwear</i> (Aspheric)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
<i>Airwear Essential Blue Series®</i> (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
Plastic (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
<i>Airwear Transitions Signature GEN 8™</i> Gray (Spherical)**	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
<i>Plastic Transitions Signature GEN 8</i> Gray (Spherical)**	1.50	+2.00 to -5.00 up to -2.00 cyl	65/70 mm

* -10.25 to -12.00 available in sphere only.

** Available May 24, 2022.

FINISHED Single Vision

CRIZAL® EASYPRO			
MATERIAL	INDEX	POWER RANGE	SIZE
<i>Thin & Lite 1.67 (Aspheric)</i>	1.67	+6.00 to -12.00 up to -2.00 cyl	65/70/75/80 mm
<i>Airwear (Spherical)</i>	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
<i>Airwear Transitions Signature GEN 8 Gray (Spherical)</i>	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
<i>Plastic (Spherical)</i>	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
<i>Plastic Transitions Signature GEN 8 Gray (Spherical)</i>	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70 mm

** -10.25 to -12.00 available in sphere only.

FINISHED Single Vision**CRIZAL EASY[®]**

MATERIAL	INDEX	POWER RANGE	SIZE
Airwear [®] (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
Plastic (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm

CRIZAL[®] KIDS UV™

MATERIAL	INDEX	POWER RANGE	SIZE
Airwear Jr. (Spherical)	1.59	+6.00 to -6.00 up to -2.00 cyl	60 mm


* -10.25 to -12.00 available in sphere only.

† Crizal Kids UV lenses are bundled products. They come as a single vision polycarbonate lens, coated with the respective Crizal[®] No-Glare lens treatment.

INTRODUCING

NEW Crizal® Sapphire™ HR

THE BEST ANTI-REFLECTIVE COATING OVERALL⁽¹⁾



High Surface Density
Process™

SUPER HYDROPHOBIC

For smudge resistance and a superior ability to repel dirt and water.

ANTI-STATIC

Embedded in the AR multi-layers to repel dust and keep lenses cleaner longer.



360° Multi-angular
Technology™

AN EMBEDDED NANOLAYER

A multi-angular (α) criterion that quantifies the overall anti-reflective efficiency (0° - 45°) to reduce reflections whatever the light direction (front and back sides of the lens)



High-Resistance
Technology™

NEW INNOVATIVE RECIPE

A combination of four specific oxides known for their resistance against dust, scratches and temperature.



E-SPF® INDEX

Crizal coatings offer the most comprehensive daily protection against the invisible and often irreversible dangers of UV reflection in a clear lens.⁽³⁾

* Stack order for demonstration only

(1) External laboratory tests and internal technical tests - 2020. Compared to the competitor most known lens-brands by consumers (2019 external brand tracking in 11 countries). The word 'overall' refers to anti-reflective coatings important criteria, ranked through an external quantitative consumer study - 2019.

(2) Number of layers for illustration purpose.

(3) E-SPF® is an index rating the overall UV protection of a lens. E-SPF® was developed by Essilor International and endorsed by third party experts. The E-SPF® index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear). *Except for Orma®: E-SPF® 10 index.

Crizal® | AN INVISIBLE SHIELD OF PROTECTION*

CORE RANGE						
	REDUCTION OF REFLECTIONS	BACKSIDE UV PROTECTION 	EASINESS TO CLEAN 	DUST REPELLENCE Anti-Static Layer	DURABLE RESISTANCE 	OPTIMAL TRANSPARENCY  Multi-angular Technology
NEW Crizal® Sapphire™ HR						
Crizal® Rock™						
NEW Crizal® Easy™ Pro						
NEW Crizal® Easy®						
	REDUCTION OF REFLECTIONS	BACKSIDE UV PROTECTION	EASINESS TO CLEAN	DUST REPELLENCE	HARMFUL BLUE LIGHT ¹ FILTRATION	
Crizal® Prevencia®						
Crizal® Kids UV	Comfort designed for kids				Comfort for kids + reduces exposure to Harmful Blue Light ¹	

* Protection defined by resistance to glare, scratching, smudging, and comprehensive UV protection on all lens materials except 1.50 clear plastic. (Internal R&D measurements - 2020 - Except clear 1.50 index)

** E-SPF® is an index rating the overall UV protection of a lens. E-SPF was developed by Essilor International and endorsed by third party experts. The E-SPF index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear).

(1) Harmful Blue Light is the blue-violet wavelengths found between 415–455nm on the light spectrum believed most toxic to retinal cells. Harmful Blue Light protection only available with Crizal Prevencia No-Glare lenses. Crizal Prevencia No-Glare lenses filter at least 20% of Harmful Blue Light.

CRIZAL® NO-GLARE LENSES

Worn and trusted by millions of people worldwide.

For the most current availability chart, visit **ESSILORUSA.COM**



©2022 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of their respective owners. *Transitions*, *Transitions Signature*, *XTRActive* are registered trademarks, and the *Transitions* logo are trademarks of Transitions Optical, Inc. used under license by Transitions Optical Ltd. GEN 8 is a trademark of Transitions Optical Ltd. Photochromic performance is influenced by temperature, UV exposure and lens material. These products may be protected by one or more patents listed at www.essilorusa.com/patents.
264421_PRO_ZAL SHK/ECST 10/22

