






Crizal®

Crizal®

FOR THE CLEAREST VISION POSSIBLE™

CRIZAL® NO-GLARE LENSES
PRODUCT GUIDE

With Crizal No-Glare lenses, patients also get a two-year warranty.*

		Crizal [®] SAPPHIRE [®] 360° UV	Crizal [®] PREVENIA [®]	Crizal [®] AVANCE [®] UV	Crizal [®] ALIZE [®] UV	Crizal [®] easy UV
CLARITY	GLARE REDUCTION Virtually reflection-free for better vision during the day and safer driving at night.	BEST-IN-CLASS CLARITY	EXCEPTIONAL CLARITY			
		VERY FAINT BLUE HUE	BLUE-VIOLET REFLECTION	FAINT GREEN HUE		
DURABILITY	DUST REPELLENT Prevents particles from adhering to the lens for less frequent cleaning. WATER REPELLENT Repels moisture for convenient outdoor wear. SCRATCH AND SMUDGE RESISTANCE Protects lenses from the wear and tear of everyday use and cleaning. Reduces smears and fingerprints for clarity without distractions.	ANTI-PARTICULATE LAYER FOR PROTECTION AGAINST DUST AND DIRT				
		SUPERIOR CLEANABILITY: HIGH SURFACE DENSITY™ TECHNOLOGY MAKES LENS SUPER HYDROPHOBIC + THE PATENTED PAD-CONTROL SYSTEM™				EASY TO CLEAN WITH PATENTED PAD-CONTROL SYSTEM™
		CONTACT ANGLE OF 116°				CONTACT ANGLE OF 112°
		2X MORE SCRATCH RESISTANCE WITH SR BOOSTER™ LAYER†		BETTER SCRATCH RESISTANCE THAN ORDINARY SCRATCH-COATED LENSES		
PROTECTION	UV PROTECTION Essential daily protection against UV light. HARMFUL BLUE LIGHT** PROTECTION Crizal® Prevenia® selectively filters out Harmful Blue Light for constant protection.					
		100% UV PROTECTION ON FRONT SIDE OF LENS*** WITH ADDITIONAL PROTECTION FROM UV REFLECTION OFF BACK SIDE OF LENS				
			FILTRATION OF UP TO 20% OF HARMFUL BLUE-VIOLET LIGHT			



Transitions™

Crizal®

VARILUX®

Eyezen™

Xperio[®] UV

*Lenses are guaranteed for up to two years from original order date from product failure resulting in atypical scratching or crazing from normal use.

Harmful Blue Light is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells. *With the exception of clear plastic (1.50) substrates. †As compared to Crizal Alize.
©2019 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries.
These products may be protected by one or more patents listed at www.essilorusa.com/patents. 54930_PRO_ZAL SHK/SSX 10/19

With Crizal lenses

Without Crizal lenses

Simulated Images

Crizal lenses combine multiple layers and technologies to offer the best combination of clarity, durability, and protection.

FOR THE CLEAREST VISION POSSIBLE™

CRIZAL FOR KIDS



Crizal *SAPPHIRE* 360° UV

REDUCES GLARE
RESISTS SCRATCHES
REPELS DUST
REPELS WATER
RESISTS SMUDGES
OPTIMAL UV PROTECTION
Best in Clarity, Scratch Resistance, Cleanability, and UV Protection

Crizal *PREVENCIA*

REDUCES GLARE
RESISTS SCRATCHES
REPELS DUST
REPELS WATER
RESISTS SMUDGES
HARMFUL BLUE LIGHT* FILTRATION
OPTIMAL UV PROTECTION
Best in UV Protection, Harmful Blue Light Filtration, and Cleanability

Crizal *AVANCE* UV™

REDUCES GLARE
RESISTS SCRATCHES
REPELS DUST
REPELS WATER
RESISTS SMUDGES
OPTIMAL UV PROTECTION
Best in Cleanability and Scratch Resistance

Crizal *ALIZE* UV™

REDUCES GLARE
RESISTS SCRATCHES
REPELS WATER
RESISTS SMUDGES
COMPREHENSIVE UV PROTECTION
Best in Cleanability

Crizal *easy* UV™

REDUCES GLARE
RESISTS SCRATCHES
REPELS WATER
RESISTS SMUDGES
COMPREHENSIVE UV PROTECTION
Everyday Essential No-Glare Lens

Crizal for Kids No-Glare lenses are a no-glare treatment, polycarbonate lens material package. They are **impact resistant** and **virtually unbreakable**—perfect for kids’ active lives!

Crizal *PREVENCIA* kids

REDUCES GLARE
RESISTS SCRATCHES
REPELS WATER
RESISTS SMUDGES
HARMFUL BLUE LIGHT* FILTRATION
OPTIMAL UV PROTECTION
Best in UV Protection and Harmful Blue Light Filtration for Kids

Crizal *kids* UV™

RESISTS SCRATCHES
REDUCES GLARE
REPELS WATER
RESISTS SMUDGES
COMPREHENSIVE UV PROTECTION
Everyday Essential No-Glare Lens for Kids

Image represents the Crizal Sapphire® 360° UV No-Glare lens treatment stack. Thickness and order of layers is for demonstration purposes only.

In a clear lens, E-SPF® is an index rating of the overall UV protection of a lens. E-SPF® was developed by Essilor International and endorsed by third-party experts. The E-SPF® index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear). *Harmful Blue Light is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

Recommended for ages 5–12. Crizal Kids UV™ is available in polycarbonate lens material only.