

VARILUX®  
#1

PROGRESSIVE LENS BRAND  
RECOMMENDED BY OPTOMETRISTS  
AND OPTICIANS\*

# VARILUX®

## LENSES AVAILABILITY

for Eyecare Professionals

Effective APRIL 2022



\* Survey conducted in 2018 by an independent third party. Sponsored by Essilor. Results were reported by independent opticians and optometrists.

VARILUX® X DESIGN		VARILUX® X FIT											VARILUX® X 4D			
		LENS ID	LENS ID											LENS ID		
														4D		
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00	
THIN & LITE®	★ Clear	1.67											-12.00 to +9.00	0.75 to 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-12.00 to +9.00	0.75 to 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00	
	★ Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00	
	Transitions® Signature® GEN 8™**	1.74	✓	✓									-14.00 to +13.00	0.75 to 4.00D	-6.00	

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors


✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.


**VARILUX® PHYSIO® W3+**



**VARILUX® PHYSIO® W3+ FIT**



**VARILUX® PHYSIO® W3+ EYECODE™**



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓		✓							-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™**	1.74	✓										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

**MFH=14mm**

**NOTE:** \*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**VARILUX® COMFORT MAX**



**VARILUX® COMFORT MAX FIT**



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive® new generation	1.53	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive® new generation	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
Transitions® XTRActive® new generation	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00	

**MFH=14mm**

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

**SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.**

**VARILUX PHYSIO DRx™** LENS ID  
H → e

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00	
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.50	✓											-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.59	✓											-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00	
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.67	✓											-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74												-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
Transitions® Signature® GEN 8™**	1.74	✓											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00	

**MFH=17mm**

**NOTES:** All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

**VARILUX PHYSIO SHORT DRx™** LENS ID  
→ e  
S

RX AVAILABILITY FOR VARILUX PHYSIO SHORT DRX MATCHES VARILUX PHYSIO DRX, EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

**MFH=14mm**

VARILUX COMFORT DRx™ 

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.50	✓											-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓	SM	SM	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.59	✓											-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓	SM	SM	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.67	✓											-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	★ Clear**	1.74												-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00
Transitions® Signature® GEN 8™**	1.74	✓											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

\*ADD Power in steps of 0.25D.

\*\*THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

\*\*\*For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

VARILUX COMFORT SHORT DRx™ 

RX AVAILABILITY FOR VARILUX COMFORT SHORT DRX MATCHES VARILUX COMFORT DRX, EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

MFH=14mm

