



CRIZAL[®] Rx AND STOCK LENS AVAILABILITY CHART

MAY 2018

TABLE OF CONTENTS

Crizal® Rx Lens Availability

Varilux® progressive lenses.....	3
DEFINITY® and Essilor PALs progressive lenses	4
Eyezen™+, Essilor Single Vision, Bifocals, Trifocals Mid Distance/Visual Fatigue, <i>Crizal</i> for Kids lenses	5
Shamir Insight, Inc. lenses.....	6
Shore lenses	7
Signet Armorlite lenses	8
KBco, ZEISS, SOLA, AO and Seiko lenses	9
Younger Optics, X-Cel lenses and Specialty Items	10
<i>Crizal</i> stock lenses	11, 12
<i>Crizal</i> lens benefits	13
Crizal Sapphire® 360° UV lens 9-step process.....	14
<i>Crizal</i> product family.....	15



Crizal® Live Life in the Clear®		PLASTIC 1.50					TREXA®				AIRWEAR®					THIN & LITE®													
		Clear ¹	Essilor Polarized ^{**4}	Transitions ^{*4}			Clear	Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage	1.60				1.67 ²				1.74 ³				
				Signature VII	XTRActive®	Vantage®											Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	
PROGRESSIVE LENSES																													
VARILUX® DIGITAL*	Varilux® X 4D	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® X Fit	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® X Design	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® S 4D	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® S Fit	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® S Design	■	●	●	●	■	■	●	■	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® Physio® W3+ eyecode™	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® Physio® W3+ Fit	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux® Physio® W3+	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	Varilux Comfort® W2+ eyecode™	■	■	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	Varilux Comfort® W2+ Fit	■	■	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	Varilux Comfort® W2+	■	■	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	Varilux Physio DRx™	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	Varilux Physio Short DRx™	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Varilux Comfort DRx™	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
Varilux Comfort Short DRx™	■	●	●	●	■	■	●	●	■	■	■	●	●	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
VARILUX®	Varilux Ipseo IV eyecode™	●										●																	
	Varilux Ipseo IV™	●										●																	
	Varilux Physio Enhanced eyecode™	●	●	●		●		●	●			●	●	●	●		●		●				●		●				
	Varilux Physio Enhanced Fit™	●	●	●		●		●	●			●	●	●	●		●		●				●		●				
	Varilux Physio Enhanced Azio™	●	●			●		●	●			●	●	●	●		●		●				●		●				
	Varilux Physio Enhanced India™	●	●			●		●	●			●	●	●	●		●		●				●		●				
	Varilux Physio Enhanced™	●	●	●		●		●	●	●		●	●	●	●		●		●				●		●				
	Varilux Physio 360°®	●	●									●	●										●						
	Varilux Physio Short 360°™											●											●						
	Varilux® Physio®	■	●	●		■		●				■	●	●	●		●		●				●						
	Varilux Physio Short™	■		●		■						■		●			●		■		●								
	Varilux Comfort Enhanced Fit™	●		●	●							●	●	●	●		●		●		●								
	Varilux Comfort Enhanced™	●		●	●							●	●	●	●		●		●		●								
	Varilux Comfort®	■		●	●							■	●	●	●		●		■		●								
	Varilux Comfort Short™	■		●	●							■	●	●	●		●		■		●								
	Varilux Ellipse 360°™	●										●	●						●				●						
	Varilux® Ellipse®	●										●	●						●				●						
	Varilux Sport™												●																
Varilux® Stylistic®	●	●	●	●	●						●	●	●	●	●			●	●	●	●								

1 E-SPF® of 35 with Crizal Sapphire® 360° UV. and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.
2 All 1.67 materials can only be tinted in Gray and Brown.
3 *Thin & Lite* 1.74 lenses come systematic with your choice of *Crizal® Previncia®*, *Crizal Sapphire 360° UV*, *Crizal Avancé UV* or *Crizal Alizé UV*.
4 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

■ Available in Essential Blue Series™ lenses.
● *Crizal Sapphire® 360° UV* available on SFSV Airwear Colors.
+ *Crizal Previncia* is not available on polarized lenses.
▲ *Crizal SunShield UV™* available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018

Better Sight.
Better Life.



PLASTIC 1.50				TREXA [*]			AIRWEAR [†]			THIN & LITE [‡]											
Clear ¹	Essilor Polarized ^{***▲}	Transitions ^{§4}			Clear	Essilor Polarized ^{***}	Transitions ⁴			1.60				1.67 ²				1.74 ³			
		Signature [®] VII	XTRActive [®]	Vantage [®]			Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ^{***}	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ^{***}	Signature VII	XTRActive	Vantage	Clear ³	Essilor Polarized [*]

PROGRESSIVE LENSES

DEFINITY	DEFINITY [®] 3 PLUS	•		•				•		•					•		•					
	DEFINITY [®] 3	•		•				•		•					•		•					
	DEFINITY [®]	•	•	•				•		•					•		•					
	DEFINITY SHORT [™]	•	•	•				•		•					•		•					
ESSILOR PALS	Essilor Ideal [®] Fit [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal [®] Advanced Fit [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal [®]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal Short [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal Advanced [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal Advanced Wrap [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal [®] Sport	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal [®] Driving	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Ideal [®] Computer [™]	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Accolade [®]	•		•				•	•	•	•	•	•	•	•	•	•	•	•	•		•
	Accolade Freedom [™]	•		•				•	•	•	•	•	•	•	•	•	•	•	•	•		•
	SmallFit [®]	•		•				•	•	•	•			•	•	•	•	•	•	•		•
	SmallFit [®] Digital	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Ovation [®]	•		•				•	•	•	•			•	•	•	•	•	•	•		•
	Ovation [®] Digital	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Essilor Natural [®]	•		•				•	•	•	•			•	•	•	•	•	•	•		•
	Essilor Natural [®] Digital	•	•	•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Adaptar [®]	•						•							•							•
	Adaptar [®] Digital	•		•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•
	Adaptar [®] Short Digital	•		•	•	•		■	•	•	•	•	•	•	•	•	•	•	•	•		•

- 1 E-SPF[®] of 35 with Crizal Sapphire[®] 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.
2 All 1.67 materials can only be tinted in Gray and Brown.
3 *Thin & Lite* 1.74 lenses come systematic with your choice of *Crizal[®] Prevencia[®]*, *Crizal Sapphire 360° UV*, *Crizal Avancé UV* or *Crizal Alizé UV*.
4 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

- Available in Essential Blue Series[™] lenses.
* *Crizal Sapphire[®] 360° UV* available on SFSV Airwear Colors.
+ *Crizal Prevencia* is not available on polarized lenses.
▲ *Crizal SunShield UV* available on polarized and tinted lenses.

Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018

Better Sight.
Better Life.



PLASTIC 1.50				TREXA [*]			MID	AIRWEAR [†]			THIN & LITE [‡]											
Clear ¹	Essilor Polarized ^{**▲}			Clear	Essilor Polarized ^{**}			Clear	Essilor Polarized ^{**}			1.60			1.67 ²			1.74 ³				
	Transitions ⁴				Transitions ⁴				Transitions ⁴			1.60			1.67 ²			1.74 ³				
	Signature [®] VII	XTRActive [®]	Vantage [®]		Signature VII	XTRActive [®]	Vantage [®]		Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage		
																					Transitions ⁴	
Essilor Polarized ^{**▲}			Essilor Polarized ^{**}			Essilor Polarized ^{**}			1.60			1.67 ²			1.74 ³							
Signature [®] VII			Signature VII			Signature VII			Signature VII			Signature VII			Signature VII			Signature VII				
XTRActive [®]			XTRActive [®]			XTRActive [®]			XTRActive			XTRActive			XTRActive			XTRActive				
Vantage [®]			Vantage [®]			Vantage [®]			Vantage			Vantage			Vantage			Vantage				



OTHER LENSES

SINGLE VISION	Eyezen [™] + Single Vision	•	•	•	•	•		•		•				•	•	•	•	•		•		•		
	Crizal Kids UV [™]							•						•										
	Crizal [®] Prevencia [®] Kids							•						•										
	SV	•	•	•	•	•		■	•	•	•	•		•	■	•	•	•	•	•	•	•	•	•
	SV AS	•		•				■		•				■	•	•	•	•		•	•	•	•	•
	Essilor 360 [°] ™							•		•	•			•						•	•	•	•	•
	Essilor eyecode [™] Single Vision							•		•	•			•						•	•	•	•	•
	Essilor FIT [™] Single Vision							•		•	•			•						•	•	•	•	•
	Essilor Azio [™] Single Vision							•		•				•						•	•	•	•	•
	Essilor Stylistic [™] WRAP Single Vision	•	•	•	•	•								•	•	•	•	•	•					
BIFOCALS	Essilor Sport Solutions [™] 5													•										
	FT28	•	•	•	•			•		•				•	•	•	•			•	•			
	FT28 AS																			•				
	FT28 Essilor Bifocal AB [™]	•		•										•						•				
	FT35	•	•	•	•			•		•				•	•	•				•	•			
	FT45	•																						
TRIFO-CALS	ROUND 22	•																						
	7 x 28	•	•	•				•	•		•			•						•	•			
MID DISTANCE	8 x 35	•	•	•				•		•				•						•	•			
	Essilor Anti-Fatigue [™]	•												•						•				
	Essilor Computer [™]													•										
MID DISTANCE	Interview [®]	•																						

- 1 E-SPF[®] of 35 with Crizal Sapphire[®] 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.
2 All 1.67 materials can only be tinted in Gray and Brown.
3 *Thin & Lite* 1.74 lenses come systematic with your choice of *Crizal[®] Prevencia[®]*, *Crizal Sapphire 360° UV*, *Crizal Avancé UV* or *Crizal Alizé UV*.
4 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.
5 *Essilor Sport Solutions* has 6 different lens colors available in tint and polarized.

- Available in Essential Blue Series[™] lenses.
* *Crizal Sapphire 360° UV* available on SFSV Airwear Colors.
+ *Crizal Prevencia* is not available on polarized lenses.
▲ *Crizal SunShield UV* available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



PLASTIC 1.50				TRIVEX			POLY			HI-INDEX										Younger Drivewear [™] 1.50 / Polycarbonate
Clear ¹	Essilor Polarized ^{++▲}	Transitions ^{*4}			Clear	Essilor Polarized ⁺⁺	Transitions ⁴			1.60			1.67 ²			1.74 ³				
		Signature ^{* VII}	XTRActive [*]	Vantage [*]			Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	Clear ²	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	
Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage						

OTHER LENSES																								
SHAMIR	Shamir SFSV	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph III® Fixed 11, 13, 15, 18 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph III Variable	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Attitude III Fashion 15 mm, 18 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Attitude III® Sport 18 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph III Single Vision	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Attitude III Single Vision	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II® Fixed 11, 13, 15, 18 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II Variable	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II Attitude 18 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II Attitude Short 15 mm	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II Office	●									●					●								
	Autograph II Single Vision	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Autograph II Single Vision Attitude™	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●	●
	Element™ 19 mm	●	●	●			■		●	●		●	●	●	●		●		●	●	●	●		
	Element Short 16 mm	●	●	●			■		●	●		●	●	●	●		●		●	●	●	●		
	Spectrum™ SV	●	●	●	●	●	■	●	●	●		●	●	●	●	●	●	●	●	■	●	●	●	●
	Spectrum 14, 16, 18 mm	●	●	●	●	●	■	●	●	●		●	●	●	●	●	●	●	●	■	●	●	●	●
	Spectrum+ 14, 16, 18 mm	●	●	●	●	●	■	●	●	●		●	●	●	●	●	●	●	●	■	●	●	●	●
	Relax™	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●	●	●	●	■	●	●	●	●
	FirstPAL™	●	●	●	●	●						●	●	●	●	●	●	●	●	●	●			●
	Golf™						■	●	●	●	●	●	●	●	●	●	●	●	●					●
	Golf - Single Vision™						■	●	●	●	●	●	●	●	●	●	●	●	●					●
	InTouch™	●	●	●	●	●	■	●	●	●		●	●	●	●	●	●	●	●	■	●	●	●	●
	WorkSpace™	●	●	●	●		■		●	●		●		●		●	●		●	■				
	Computer™	●	●	●	●		■		●	●		●		●		●	●		●	■				
	Shamir Duo™	●	●	●	●	●	■	●	●	●	●	●	●	●	●	●		●	●	■		●	●	●

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.

2 All 1.67 materials can only be tinted in Gray and Brown.

3 *Thin & Lite* 1.74 lenses come systematic with your choice of *Crizal® Previncia®*, *Crizal Sapphire 360° UV*, *Crizal Avancé UV* or *Crizal Alizé UV*.

4 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

■ Available in Blue Zero™.

* *Crizal Sapphire® 360° UV* available on SFSV Airwear Colors.

+ *Crizal Previncia* is not available on polarized lenses.

▲ *Crizal SunShield UV* available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



PLASTIC 1.50				TRIVEX ¹			POLY			HI-INDEX									
Clear ¹	Essilor Polarized ^{**+▲}	Transitions ⁻³			Clear	Essilor Polarized ^{**}	Transitions ⁻³			1.60				1.67 ²					
		Signature [®] VII	XTRActive [*]	Vantage [*]			Signature [®] VII	XTRActive [*]	Vantage [*]	Clear	Essilor Polarized ^{**}	Transitions ⁻³			Clear ²	Essilor Polarized ^{**}	Transitions ⁻³		
												Signature [®] VII	XTRActive [*]	Vantage [*]			Signature [®] VII	XTRActive [*]	Vantage [*]

OTHER LENSES																											
SHORE	Shore SV Spherical	•	•				•					•	•				•					•	•				
	Shore SFSV AS											•										•					
	Shoreview® ES	•		•								•	•	•			•					•					
	Shoreview Mini ES	•		•								•	•	•			•					•					
	ShoreExec																										
	Shoreview	•					•					•	•				•										
	Shoreview Mini	•					•					•					•										
	Shoreview Digital	•		•						•			•	•	•			•			•		•	•	•		
	Shoreview Digital Mini	•		•						•			•	•	•			•			•		•	•	•		
	Shoreview Digital Advanced	•	•	•	•	•	•			•		•	•	•	•	•	•	•		•		•	•	•			
	Shoreview Digital Advanced Mini	•	•	•	•	•	•			•		•	•	•	•	•	•	•		•		•	•	•			

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.

2 All 1.67 materials can only be tinted in Gray and Brown.

3 *Transitions*[®] products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

• *Crizal Sapphire® 360° UV* available on SFSV Airwear Colors.

+ *Crizal® Previncia*[®] is not available on polarized lenses.

▲ *Crizal SunShield UV*[®] available on polarized and tinted lenses.



PLASTIC 1.50					TRIVEX®			POLY				HI-INDEX						Younger Drivewear™ 1.50 / Polycarbonate
Clear¹	Essilor Polarized**⁴	PhotoViews™	Transitions⁴		Clear	Essilor Polarized**	Photo Views™	Transitions⁴		1.60			1.67²			1.74³		
			Signature® VII	XTRActive®				Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Clear²	Essilor Polarized**	Signature VII	XTRActive	
	Transitions⁴		Clear	Essilor Polarized**	Photo Views™	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized**	Signature VII	XTRActive	Vantage	Clear³	Essilor Polarized			
	Transitions⁴																	
	Transitions⁴																	

OTHER LENSES

[illegible]

1 E-SPF* of 35 with Crizal Sapphire® 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.

2 All 1.67 materials can only be tinted in Gray and Brown.

3 *Thin & Lite 1.74* lenses come systematic with your choice of *Crizal® Prevencia®*, *Crizal Sapphire 360° UV*, *Crizal Avancé UV* or *Crizal Alizé UV*.

4 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

* *Crizal Sapphire® 360° UV* available on SFSV Airwear Colors.

+*Crizal Previncia* is not available on polarized lenses.

▲ *Crizal SunShield UV*® available on polarized and tinted lenses.



PLASTIC 1.50				TRIVEX®			MID	POLY			HI-INDEX									
Clear¹	Essilor Polarized**⁴	Transitions³			Clear	Essilor Polarized**	1.55 Essilor Polarized**	Clear	Essilor Polarized**	Transitions³			1.60			1.67²				
		Signature VII	XTRActive*	Vantage*						Signature VII	XTRActive	Vantage	Signature VII	XTRActive	Vantage	Clear²	Essilor Polarized**	Transitions³		
																		Signature VII	XTRActive	Vantage

OTHER LENSES

KBco	SV		•					•							•						•				
	SV AS																				•				
	EOS™ PAL		•																			•			
	iRx® PRO		•						•													•			
	iRx® SHORT								•																
	Xplore™		•									•													
	FT28		•							•		•									•	•			
	FT35		•							•		•									•	•			
	7 x 28		•							•											•	•			
	8 x 35		•							•											•	•			
ZEISS	ZEISS GT2®	•		•								•								•		•			
	ZEISS GT2® Short	•		•								•								•		•			
	ZEISS Gradal Top	•		•																					
SOLA	SOLA Compact ULTRA™ HD	•		•																•		•			
	SOLA MAX®	•		•																					
	SOLAOne® HD	•		•																•		•			
AO	AO® Compact	•		•																					
SEIKO	SFSV																			•					
	Proceed II																			•					
	Proceed III																			•					

1 E-SPF® of 35 with Crizal Sapphire® 360° UV and E-SPF of 25 with all other *Crizal* products made with any lens material other than clear 1.50 plastic.

2 All 1.67 materials can only be tinted in Gray and Brown.

3 *Transitions* products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

*Crizal Sapphire® 360° UV available on SFSV Airwear Colors.

+ *Crizal Prevencia* is not available on polarized lenses.

▲ *Crizal SunShield UV*® available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



PLASTIC 1.50				TRIVEX [*]			MID	POLY			HI-INDEX									Younger Drivewear [®] 1.50 / Polycarbonate						
Clear ¹	Essilor Polarized ^{++▲}	Transitions ⁻³		Clear	Essilor Polarized ⁺⁺	Transitions ⁻³		1.55 Essilor Polarized ⁺⁺	Clear	Essilor Polarized ⁺⁺	Transitions ⁻³		1.60			1.67			1.74 ²							
		Signature [®] VII	XTRActive [®]			Vantage [®]	Signature VII				XTRActive	Vantage	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ⁺⁺	Transitions ⁻³		Signature VII	XTRActive	Vantage	Clear ³	Essilor Polarized [®]	

OTHER LENSES																								
YOUNGER	SV Spherical						•																	•
	SV Aspherical						•																	
	FT28												•											
	Image [®]	•	•	•			•																	•
X-CEL	SFSV	•		•	•	•	•		•	•														
	Asp-SFSV	•			•	•	•																	
	D-25	•																						
	D-28	•		•	•		•		•															
	D-35	•		•			•		•															
	7X25	•																						
	7X28	•		•	•		•		•															
	7X35	•																						
	8X35	•		•			•		•															
	RD-22	•																						
	RD-24			•			•		•															
	Blended Rd	•																						
	Freedom ID Pal			•			•		•	•	•													
SPECIALTY ITEMS (on UNCOATED LENSES ONLY)																								
	Slaboff, High Adds, Bi-Concave, etc...	•																						

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360° UV and E-SPF of 25 with all other Crizal products made with any lens material other than clear 1.50 plastic.

2 Thin & Lite 1.74 lenses come systematic with your choice of Crizal[®] Previncia[®], Crizal Sapphire 360° UV, Crizal Avancé UV or Crizal Alizé UV.

3 Transitions products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

• Crizal Sapphire[®] 360° UV, not available on polarized or tinted lenses.

+ Crizal Previncia is not available on polarized lenses.

▲ Crizal SunShield UV[®] available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



FINISHED SINGLE VISION

CRIZAL [®] PREVENCIA [®]			
MATERIAL	INDEX	POWER RANGE	SIZE
Airwear [®] (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl Available to a -6.00	
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
		-2.00 cyl Available to a -4.00	

CRIZAL AVANCÉ UV [™]			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite [®] 1.74 (Aspheric)	1.74	-3.00 to -12.00 up to -2.00 cyl*	70/75/80 mm
		-2.00 cyl available to a -10.00	
Thin & Lite 1.67 (Aspheric)	1.67	+6.00 to -10.00 up to -2.00 cyl	65/70/75/80 mm
		2.00 cyl available to a -10.00	
Thin & Lite 1.67 Transitions [®] Signature [®] VII Gray (Aspheric)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75/80 mm
		-2.00 cyl available to a -8.00	
Thin & Lite 1.60 (Aspheric)	1.60	+6.00 to -10.00 up to -2.00 cyl	65/70/75 mm
		-2.00 cyl available to a -10.00	
Airwear [®] (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -6.00	
Airwear Transitions Signature VII Gray & Graphite Green (Spherical)	1.59	4.00 to -6.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -4.00	
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
		-2.00 cyl available to a -4.00	
1.50 Transitions Signature VII Gray (Spherical)	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -2.00	

* -10.25 to - 12.00 available in sphere only.

FINISHED SINGLE VISION

CRIZAL ALIZÉ UV™			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite® 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75 mm
		-2.00 cyl available to a -6.00	
Thin & Lite 1.67 (Aspheric)	1.67	+6.00 to -12.00 up to -2.00 cyl	65/70/75/80 mm
		-2.00 cyl available to a -10.00*	
Thin & Lite 1.60 (Spherical)	1.60	+6.00 to -10.00 up to -2.00 cyl	65/70/75 mm
		-2.00 cyl available to a -10.00	
Airwear® (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -6.00	
Airwear® Transitions® Signature® VII Gray (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -4.00	
TREXA® (Spherical)	1.53	+4.00 to -6.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -6.00	
Plastic Flat (Spherical) ¹	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
		-2.00 cyl available to a -4.00	
1.50 Transitions® Signature® VII Gray (Spherical)	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -2.00	
CRIZAL EASY UV™			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75 mm
		-2.00 cyl available to a -6.00	
Thin & Lite 1.67 (Aspheric)	1.67	-2.00 to -12.00 up to -2.00 cyl	70/75/80 mm
		-2.00 cyl available to a -10.00*	
Airwear (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70 mm
		-2.00 cyl available to a -6.00	
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl	70 mm
		-2.00 cyl available to a -4.00	

* -10.25 to -12.00 available in sphere only.

CRIZAL® NO-GLARE LENSES PROVIDE THE CLEAREST VISION POSSIBLE



Reduce glare

Glare, poor night vision, and reflection are major concerns for patients

- Less reflection during the day and safer nighttime driving



Resist smudges

Lenses that smudge easily are a major source of frustration

- Fewer smears and fingerprints mean less visual distraction



Resist scratches

Scratch resistance is the benefit patients want most in anti-reflective lenses

- Protection against everyday wear and tear, with no scratches from cleaning



Repel dust

Patients are bothered by dust, dirt, and lashes clinging to static on the inside of lenses

- *Crizal* repels these particles, so lenses stay cleaner longer



Repel water

Patients want lenses that stay clean in outdoor environments


- More convenient to wear outdoors



Comprehensive UV protection

UV protection is needed every day, but many patients don't know this

CRIZAL SAPPHIRE® 360° UV LENS 9-STEP PROCESS



E-SPF is an index rating the overall UV protection of a lens. E-SPF was developed by Essilor International and endorsed by third party experts. The E-SPF index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear).

Image represents *Crizal Sapphire 360° UV* stack
*Thickness of layers for demonstration purposes only. The image is not representative of the actual thickness. Total Crizal lens process thickness is less than 0.5p. Crizal.com 14 of 16

CRIZAL® PRODUCT FAMILY



- Crizal easy UV™**
 - More durable than ordinary lenses
 - Easier to clean than ordinary lenses
- Crizal ALIZÉ UV™**
 - All the benefits of *Crizal Easy UV* lenses, PLUS:
 - Superior protection against smudges, oil, and water
 - Lenses stay cleaner longer
- Crizal AVANCÉ UV™**
 - All the benefits of *Crizal Alizé UV* lenses, PLUS:
 - 2x more scratch resistance to protect lens investment
 - Repels dust and dirt for clearer vision and less cleaning
- Crizal SAPPHIRE® 360° UV**
 - All the benefits of *Crizal Avancé UV* lenses, PLUS:
 - Best-in-class transparency
 - Reduction of reflections from all around the wearer
 - Plus, UV protection has been improved with E-SPF® 35
- Crizal PREVENCIA®**
 - Selectively filters at least 20% of Harmful Blue Light*

*Harmful Blue Light is the blue-violet wavelength between 415-455nm on the light spectrum believed most toxic to retinal cells.



Applies to the Crizal lens process

CRIZAL[®] NO-GLARE LENSES

Worn and trusted by millions of people worldwide.



Optimal UV protection in a clear lens

The higher the E-SPF[®] Index, the greater the level of UV protection provided by the lens

E-SPF is an index rating the overall UV protection of a lens. *E-SPF* was developed by Essilor International and endorsed by third party experts. The *E-SPF* index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear).

For the most current availability chart, visit **Crizal.com**



Better Sight.
Better Life.

Transitions[™]

Crizal[®]

VARILUX[®]

Eyezen[™]

Xperio[™]

"©2018 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all product names, trademarks, and company names are the property of their respective owners. Essilor, Crizal, E-SPF, and the E-SPF design are trademarks of Essilor International. TRIVEX is a registered trademark of PPG Industries Ohio, Inc. Transitions is a registered trademark and the "Bisected O" design are trademarks of Transitions Optical, Inc., used under license by Transitions Optical Ltd. **Photochromic and polarization performance is influenced by temperature, UV exposure, and lens material.** These products may be protected by one or more patents listed at www.essilorusa.com/patents. 20255_PRO_ZAL SHK/ECSL 4/18"